

VOGEL GmbH & Co. KG was founded in 1952 by Mr. Wilhelm Vogel and was taken over in the second generation by Marlies and Thomas Knauber in 1980. Today, the company is still a 100 % family owned business. VOGEL produces, develops and sells complete systems of high quality, equipments and single parts. The product line includes exhaust systems for histolabs, sophisticated stainless steel constructions (for example grossing-, cutting table and working table for staining) as well as lab-equipments, storage cabinets for formalin baskets and filing cabinets as a complete system.

- Made in Germany
- Customized stainless steel
- Individual solution for your laboratory
- Laboratory equipment
- Storage systems
- Cassettes and slide printers

VOGEL
medical technology and electronic

VOGEL GmbH & Co. KG
Gottlieb-Daimler-Straße 2
D-35463 Fernwald Annerod

Fon: +49 (641) 93 231-0
Fax: +49 (641) 39 22 1

zentrale@vogel-giessen.de
www.vogel-giessen.de






EZ-Writer-6001

manual printer for cassettes or slides





Innovative and effective the EZ-Writer by VOGEL

Compact slide or cassette printer for on-demand printing. The unit offers flexibility, consistency and reliability for slide or cassette printing needs. Can be used with computer, bar code scanner, LAN-access and USB Windows-Driver. The customer has the possibility to use the printer as an online job printer or as a stand-alone unit.

Slide printer

4 lines, default font, without bar code:

4 x 16 = 64 characters

4 lines, default font, with bar code:

2 x 11 = 22 characters (with bar code)

2 x 16 = 32 characters (without bar code)

2 lines, uppercase, without bar code:

2 x 13 = 26 characters

Prints on color slides

Cassette printer

2 lines, default font, without bar code:

2 x 16 = 32 characters

2 lines, default font, with bar code:

2 x 11 = 22 characters

1 line, uppercase, without bar code:

1 x 13 = 13 characters

Accepts standard cassettes

Similarity

- On-demand printing with or without Matrix Code
- Accepts all standard size microscope color slides or cassettes
- Number key PAD for quick number entry available
- Special characters available
- Recall saved screens from memory
- Stored computer information can be downloaded
- Upper and lower case letters
- Different printing mode
- Capable of printing special characters
- Stores up to 100 different programmable formats
- Bar code scanner adaptable (Y-gate)
- Prints cassette in approximately 2–3 seconds, slides in approximately 3 seconds
- USB/LAN-access makes it easy to transfer external data from a computer/laptop to the printer
- USB Windows Driver available
- Ribbon with UV-sensitive substance and RFID-chip
- Size: 21 x 18 x 25 cm