100944_BROCHURE:Layout 1 12/10/09 6:55 PM Page 1

NEUTRA-GUARD® SAFE, EFFECTIVE, AND ECONOMICAL FORMALIN CONTROL SYSTEM



INCINERATION = ATMOSPHERIC POLLUTION

Reduce your Carbon Footprint. Why Incinerate when Eco-Friendly "Down the Drain" Disposal of your Formalin Waste is Possible with Neutralex.®

The only **Environmentally Safe** and **Certified Technology** for Neutralization of Formalin Waste.

NEUTRALEX®
Eliminates the need for off-site disposal of hazardous formalin waste and the liability associated with its regulation. Tested using EPA compliant protocol and certified for more than 10 years, **Neutralex® is the only acceptable aldehyde neutralizer that meets the State of California EPA's requirements for safe and effective drain disposal of 10% formalin waste.**

NEUTRALEX® FAMILY

4047: 16 Pre-measured powder pouches (each neutralizes 1 gal or 4 liters)

BULK

4048: 25-Kg pail (neutralizes approximately 150 liters)

4049: 9-Kg pail (neutralizes approximately 50 liters)

WHAT DOES THE CERTIFICATION MEAN? California EPA has independently tested and evaluated Neutralex®. It certifies that all waste treated with Neutralex® will not pose any hazard to public health and safety, or to the environment. **It is the only product on the market today stamped with EPA's seal of approval. If it is not certified, it is not approved for use!**

Cal/EPA Certified Technology
Less than 3 minutes of employee time
Final neutralization less than 1ppm
No pH adjustment needed
No drain clogging solids
Safe & Easy to use



NEUTRALEX®, THE ONLY CAL/EPA CERTIFIED TECHNOLOGY



NEUTRA-PADS®
Encapsulate and neutralize 10% formalin while reducing vapors during grossing, keeping your area clean and dry. Each 8 x 10" pad absorbs up to 80-ml of formalin. **New** 11 x 17" size is great for larger specimens or for lining specimen and reagent storage areas.

4295: 100 grossing pads/case (8 x 10")
4295-1: 10 grossing pads/bag (8 x 10")
4297: 15 grossing pads/bag (11 x 17")



NEUTRA-WIPES®
Use this handy wipe to rapidly absorb and neutralize small formalin spills. Quickly reduces the effects of noxious vapors. Wipes are ideal as a liner in tissue transporting trays or as a shelf liner in storage areas. Convenient pop-up dispenser.

4278: 60 wipes 12 x 12.5"



NEUTRALEX® STARTER KIT 4046
Starter Kit Includes
4 packets of Neutralex®
10-pk 8 x 10" Neutra-Pads®
60 Neutra-Wipes®
1 Reaction Container
1 Formalin Test Kit with pH Test Strips
1container Neutra-Form®

LOOK FOR THE CERTIFIED SEAL



NEUTRALEX® REACTION CONTAINER 4276
Waste collection and Formalin treatment container
Large mouth opening for easy disposal
Screw on cap for vapor containment
Easy to read measurements
9.5 L (2.5 gal) capacity



NEUTRALEX® FORMALIN TEST KIT 4050
Fast and reliable determination of pH, formaldehyde, glutaraldehyde and paraformaldehyde in wastewater. Tests with these strips are completed in seconds: dip, remove, and compare color with the color scale provided. **Each kit contains 100 formalin and 100 pH determinations.**



NEUTRA-FORM® 4304
For accidental spills, this fast, effective control agent diminishes harmful vapors and neutralizes hazardous aldehydes. One container of control agent will absorb and neutralize a large size spill of 10% formalin, or 2% glutaraldehyde. **Two containers of control agent per case.**

Scigen

333 E. Gardena Blvd
Gardena, CA 90248
(310) 324-6576
www.scigen-direct.com

Join the thousands of satisfied professionals that are using the Neutra-Guard® family of products.

Order now at www.scigen-direct.com or Call (310) 324-6576

We are CAL/EPA Certified!

© 2009 Scigen, Inc. Rev. A